



MEDICAL DEVICES & MORE

NEPHROLOGY CARE PRODUCTS



VISION

To ensure easy reach and regular availability of our products to all our customer, at all time to continue building trust among Stake holders of **SCEPTRE**.

MISSION

To serve Medical fraternity with products of international standards and create opportunities for healthy and consistent growth of our Associates, **SCEPTRE** has never lost its vision to identify the needs of Patients and serve them with most suited products.

QUALITY POLICY

By incorporating latest manufacturing techniques we at **SCEPTRE** are committed to produce best quality products with suitable prices. Fulfilment of Customer needs and exceed their expectation are our utmost priorities. This is ensured by providing environment for hassle-free work culture and uninterrupted supply of materials, components and other resources. Our goal is to further enhance markets in exports, domestic & institutional areas and improve our quality management system.



your
kind of
care

Single Lumen Hemodialysis Catheters

- **Soft Tip** reduces trauma to vessel wall upon introduction and on-site.
- **Heparin Cap** is anti-microbial & specially-designed to help avoid infection.
- **Suture Wing** are rotating type to provides secure external anchoring.
- **Lumen** is of sufficiently thick radiopaque material to ensure correct placement of catheter when visualized under x-ray.
- **Biocompatible Polyurethane** material tested for toxicity and hemocompatible i.e., it gets softer at temperature.

Product Code	Size	Standard Packing Box/Carton	Dimensions Box/Carton
301.A111EG	7.0Fr (14G) x 13.0	25/250	110x110x350mm/570x360x270mm
301.A111EJ	7.0Fr (14G) x 14.0	25/250	110x110x350mm/570x360x270mm
301.A111FG	8.0Fr (13G) x 13.0	25/250	110x110x350mm/570x360x270mm
301.A111FJ	8.0Fr (13G) x 14.0	25/250	110x110x350mm/570x360x270mm
301.A111HG	9.0Fr (12G) x 13.0	25/250	110x110x350mm/570x360x270mm
301.A111HJ	9.0Fr (12G) x 14.0	25/250	110x110x350mm/570x360x270mm

Note: Single Lumen Hemodialysis (Femoral) Catheter are also available with following components.

- 1. Dilator
- 2. Guidewire
- 3. Introducer Needle

Available in Brand 



“Catheters and other key components are developed based on research done in USA, TAIWAN & PRC”

Transducer Protector

Antibacterial hydrophobic Air-filters with both female & male luer lock. This protective hydrophobic barrier allows only sterile air passage protecting patients and equipment from cross contamination.

Coloured Housing available optionally for identification.

Custom-Printed Blister packaging available optionally.

Housing Design is robust & safe to withstand pneumatic-pulsing.



Biological Safety

2.5cm²

Material of Construction

Media : PTFE, Housing : PVC/ABS

Intel / Outlet Connctions

Female Luer Lock (FLL)

Female Luer Slip (FLS)

Male Luer Lock (MLL)

Male Luer Lock (FLL)

All Luers confirms to ISO:594/80369-7

Maximum Operation Temperature

55°C (131°F)

Pore Size

0.2 u m

Sterilization Compatibility

ETO

Biological Safety

Materials of construction pass USP Biological Reactivity Tests, in Vivo<88>

Minimum Water Breakthrough

PTFE : 1.1 bar/30 seconds

Chemical Tests Approved

PH < 1.5ml, Heavy Metal <1ppm

Potassium Permanganate - Reducing

Substances < 1.0 ml

UV Spectrum < 0.05 ABS

(based on 350nm)

Residue on Evaporation < 1 mg

Connected with Tubing and Dialysis Machine Sensor

Transducer Protector	Product Code	Standard Packing Box / Carton	Dimensions Box/Carton	Gross Wt.
HOUSING/ Semi-Rigid PVC/ABS (Medical-Grade) MEMBRANE/PTFE (GORE™ Medical Membrane) Pore Size - 0.2 u m "Sterilization Method - ETO"	302.E0000A	100/1000 50/1000	170x135x150mm/695x330x180mm 170x135x70mm/695x330x180mm	4.0 Kg

Disposable Hemodialysis Tubing Set (Pre/Post Pump)



Product Code	Standard Packing	Dimensions Carton	Gross Wt.
301.F0000A/302.F0000A/303.F0000A	24 Pcs	540x440x200mm	7.50 Kg

Large Finger Guard on injection-site.

Tubing Material is latex-free.

The Clamps (small and large) are designed with ideal ergonomic features, no additional clamps required.

Dialyser Connector is wing-shaped for better grip and is provided with a security-lock for dialyser ports.

The Caps can also be used as safety-caps on the Dialyser.

Patient Connector grip is secure and rotating-type for smooth locking to patient's vascular access.

Drip Chamber is available in all sizes and models; rigid or soft, single or multiple, inlet/outlets & conical or straight depending upon customer needs.

Recirculator Connector can connect Venous to Arterial line for re-circulation.

Can be used as attachment hook.

Can be used as flow resistance in-line during priming (closed cap vented).

Act as safety barrier against contamination at the patient-connector before connection.

No coloured parts in direct contact with blood.

Accessories such as in Infusion Line, Priming Bag & Transducer Protector are available separately or pre-connected as and if required.

Introducing

Pediatric Long Term Hemodialysis Catheters

Designed Under Technical Collaboration With Martech USA

Product Code	Size	Standard Packing Box/Carton	Dimensions Box/Carton
301.A221FM	8.0 Fr x 18.0 cm	01/24	320x235x25mm/490x330x340mm

Contents

- One** Polyurethane Radiopaque Cuffed Catheter with Injection Site Caps
- One** Nitinol Core Guide Wire With Advancer
- Two** Dilator
- One** Scalpel Blade No. 11
- One** Introducer Needle
- One** 10 ml Syringe
- Two** Heparin Locks
- One** IV Dressing
- One** Surgical Drape
- One** Silk Surgical Suture
- One** Sheath/Sheath Valved
- One** Tunneler (Trocar)
- Four** Gauze

Adult Long Term Hemodialysis Catheters

Product Code	Size (Tip to Hub)	Standard Packing Box/Carton	Dimensions Box/Carton	Gross Wt.
301.A221INS	12.5 Fr x 28.0 cm	01/24	320x235x25mm/490x330x340mm	6.25 Kg
301.A221INT	12.5 Fr x 32.0 cm	01/24	320x235x25mm/490x330x340mm	6.25 Kg
301.A221SR	14.5 Fr x 24.0 cm	01/24	320x235x25mm/490x330x340mm	6.25 Kg
301.A221SS	14.5 Fr x 28.0 cm	01/24	320x235x25mm/490x330x340mm	6.25 Kg
301.A221ST	14.5 Fr x 32.0 cm	01/24	320x235x25mm/490x330x340mm	6.25 Kg
301.A221SV	14.5 Fr x 36.0 cm	01/24	320x235x25mm/490x330x340mm	6.25 Kg
301.A221SW	14.5 Fr x 40.0 cm	01/24	320x235x25mm/490x330x340mm	6.25 Kg

	Tip to Hub	Tip to Cuff	End to End
23cm (Tip to Cuff)	24	19	36
	28	23	40
28cm (Tip to Hub)	32	27	44



Double & Triple Lumen Hemodialysis Catheters

Double Lumen Catheter Subclavin & Curved Extensions

Subclavin Extensions	Curved Extensions	Size	Standard Packing Box/Carton	Dimensions Box/Carton	Gross Wt.
301.A121DB	301.A122DB	6.5 Fr x 8-10 cm	06/24	320x230x165mm/490x330x340mm	4.0 Kg
301.A121GC	301.A122GC	8.5 Fr x 10-11 cm	06/24	320x230x165mm/490x330x340mm	4.0 Kg
301.A121JF	301.A122JF	10 Fr x 12 cm	06/24	320x230x165mm/490x330x340mm	4.0 Kg
301.A121LJ	301.A122LJ	11.5 Fr / 12 Fr x 13.5 cm	06/24	320x230x165mm/490x330x340mm	4.0 Kg
301.A121MK	301.A122MK	11.5 Fr / 12 Fr x 15 cm	06/24	320x230x165mm/490x330x340mm	4.0 Kg
301.A121LL	301.A122LL	11.5 Fr / 12 Fr x 16 cm	06/24	320x230x165mm/490x330x340mm	4.0 Kg
301.A121LP	301.A122LP	11.5 Fr x 20 cm	06/24	320x230x165mm/490x330x340mm	4.0 Kg
301.A121MP	301.A122MP	12 Fr x 20 cm	06/24	320x230x165mm/490x330x340mm	4.0 Kg

Triple Lumen Catheter Subclavin & Curved Extensions

Subclavin Extensions	Curved Extensions	Size	Standard Packing Box/Carton	Dimensions Box/Carton	Gross Wt.
301.A131KH	301.A132LH	11.5 Fr / 12 Fr x 13.5 cm	06/24	320x230x165mm/490x330x340mm	4.0 Kg
301.A131KK	301.A132MK	11.5 Fr / 12 Fr x 15 cm	06/24	320x230x165mm/490x330x340mm	4.0 Kg
301.A131LL	301.A132LL	11.5 Fr / 12 Fr x 16 cm	06/24	320x230x165mm/490x330x340mm	4.0 Kg
301.A131KP	301.A132LP	11.5 Fr / 12 Fr x 20 cm	06/24	320x230x165mm/490x330x340mm	4.0 Kg

Custom Built Kits Available in Brands **JETLITE®** **JETFLOW®** **REIN®**
COMPOUND



Soft Tip reduces trauma to vessel wall upon introduction and on-site.

Heparin Cap is anti-microbial & specially-designed to help avoid infection.

Suture Wing are rotating type to provides secure external anchoring.

Lumen is of sufficiently thick radiopaque material to ensure correct placement of catheter when visualized under x-ray.

Biocompatible Polyurethane material tested for toxicity and hemocompatible i.e., it gets softer at temperature.

Subclavin Extensions	Curved Extensions	Size	Standard Packing Box/Carton	Dimensions Box/Carton	Gross Wt.
30I.U121DB	30I.U122DB	6.5 Fr x 8-10 cm	12	380x360x310mm	4.4 Kg
30I.U121GC	30I.U122GC	8.5 Fr x 10-11 cm	12	380x360x310mm	4.4 Kg
30I.U121JF	30I.U122JF	10 Fr x 12 cm	12	380x360x310mm	4.4 Kg
30I.U121LH	30I.U122LH	11.5 Fr / 12 Fr x 13.5 cm	12	380x360x310mm	4.4 Kg
30I.U121MK	30I.U122MK	11.5 Fr / 12 Fr x 15 cm	12	380x360x310mm	4.4 Kg
30I.U121LL	30I.U122LL	11.5 Fr / 12 Fr x 16 cm	12	380x360x310mm	4.4 Kg
30I.U121LP	30I.U122LP	11.5 Fr x 20 cm	12	380x360x310mm	4.4 Kg
30I.U121MP	30I.U122MP	12 Fr x 20 cm	12	380x360x310mm	4.4 Kg


Subclavin Extensions	Curved Extensions	Size	Standard Packing Box/Carton	Dimensions Box/Carton	Gross Wt.
301.U131LH	301.U132KH	11.5 Fr / 12 Fr x 13.5 cm	12	380x360x310mm	4.4 Kg
301.U131MK	301.U132KK	11.5 Fr / 12 Fr x 15 cm	12	380x360x310mm	4.4 Kg
301.U131LL	301.U132LL	11.5 Fr / 12 Fr x 16 cm	12	380x360x310mm	4.4 Kg
301.U131LP	301.U132KM	11.5 Fr / 12 Fr x 20 cm	12	380x360x310mm	4.4 Kg

WE ALSO MANUFACTURE MRI SAFE COAXIAL HEMODIALYSIS KIT

One	Polyurethane Radiopaque Cuffed Catheter with Injection Site Caps
One	Nitinol Core Guide Wire With Advancer
Two	Dilator
One	Scalpel Blade No. 11
One	Introducer Needle
One	5 ml Hollow Syringe
Two	Heparin Locks
One	IV Dressing
One	Surgical Drape
One	Needle Holder
One	Hole Sheet
One	Swab Stick
Three	Sterile Brushes
One	Silk Surgical Suture


Introducer Needle

- Insertion forces** extremely low.
- Easy access** & reduce tissue trauma.
- Smooth Insertion** of guide wire .
- Needle Cover** for provided to avoid damage.
- Variety of sizes** with different lengths and shapes available.



Type	Description	Product Code	Standard Packing Box/Carton	Dimensions Box/Carton
Straight	17G x 60mm	301.G0011C	100/1000	170x135x150mm/695x330x180mm
	18G x 70mm	301.G0012D		
	20G x 35mm	301.G0013A		
Y-Shape	17G x 60mm	301.G0021C	100/1000	170x135x150mm/695x330x180mm
	18G x 70mm	301.G0022D		
	20G x 35mm	301.G0023A		


Guide Wire



- Accurate Measure** for catheter tracking.
- Smooth Tip** for insertion (curved and Straight tip).
- Smooth Advancer** for Dispensing.
- Variety of sizes** with different lengths and shapes available.
- Easy Insertion** and removal.
- Kink Resistance**.
- Safety Cover** provided to avoid damage.

Product Code	Guide Wire	Tip Shape Material	Product Code	Guide Wire	Size	Tip Shape Material	Standard Packing Carton
301.L0111A	0.018" X 45	Straight/ Nitinol	301.L0121A	0.018" X 45	40-50cm (400-500mm)	Jugular/ Nitinol	75
301.L0112B	0.021" X 50	Straight/ Nitinol	301.L0122B	0.021" X 50	50-60cm (500-600mm)	Jugular/ Nitinol	75
301.L0114D	0.035" X 70	Straight/ Nitinol	301.L0124D	0.035" X 70	60-70cm (600-700mm)	Jugular/ Nitinol	75

Peritoneal Catheter



Product Code	Size
301.9111AA	8Fr x 20cm.
301.9111BC	10Fr x 20cm-32cm.
301.9111CD	12Fr x 27cm-37cm.
301.9111DE	14Fr x 47cm-57cm.
301.9111EF	16Fr x 47cm-57cm.

Features



A.V. Fistula Needle

Siliconized Needle are imported from Japan for smooth penetration.

Safety Clamps are strong & secure.

Butterfly Wings are easy to handle.

Manufacturing Process completely done in class 10,000 clean room.

Connectors are safe & compatible with all standard Hemodialysis Blood Tubing Sets.

Safety Guard provided optionally for better infection control.

Needle Gauge	Needle Length	Wing Type	Tube Length	Back Eye	Wing Color	Mini Clamp	Product Code	Standard Packing Box/Carton
14G(2.0mm)	1.0" & 1.25"	Fixed & Rotating	15.0 & 25.0 cm	With back eye & Without back eye	Purple	Yes	301.D0201A	25/500
15G(1.8mm)					Off-White	Yes	301.D0202A	25/500
16G(1.6mm)					Green	Yes	301.D0203A	25/500
17G(1.4mm)					Orange	Yes	301.D0204A	25/500

Dimensions Box/Carton	Gross Wt.
420x130x90mm/540x440x200mm	5.0 Kg
420x130x90mm/540x440x200mm	5.0 Kg
420x130x90mm/540x440x200mm	5.0 Kg
420x130x90mm/540x440x200mm	5.0 Kg



www.sceptremedica.com

Clinically Tested

ISO 9001 : 2015 & 13485 Certified

FDA Registered

Sceptre Medical Devices Pvt Ltd

Office : 6-Central Market, Ashok Vihar, PH-I, Delhi-110 052, DL, India

Works : 187, Sector-38, PH-I, HSIIDC, Rai, Sonapat-131 029, HR, India

Email : info@sceptremedica.com | **Tel :** +91 130 2368187

USA : 2100 North Line Street, J101 Lansdale, Pa 19446

Email : medical@prem4biz.com | **Tel :** +1(208)-861-6488

